



CLEANING + DISINFECTION in one step!

All-in-one solution for RGP lenses

Multi-purpose solution for the most demanding users and those who need a higher level of cleanliness.

Properties

Multi-purpose solution. GP multi users can clean, disinfect, rinse, moisten and preserve their lenses with only one solution.

The combination between HPMC and PVP increases ocular comfort so that the lenses can be used for longer.

The action of poloxamer offers extraordinary cleanliness, which, together with its anti-adhesive properties, prevents the deposits of microorganisms, hence increasing the protection against infections.

Composition

Hypromellose & PVP: lubricating agents to improve the insertion of the lens onto the eye.

Sodium chloride: which offers the necessary osmolarity.

Borax/Boric: which balances the solution's pH so that it adapts to that of the eye.

Poloxamer: is the surfactant agent that cleans the lens.

Polyhexanide: antimicrobial disinfectant.

EDTA: chelating agent, which helps to eliminate calcium deposits.